



CENTEXBEL

Rue du Travail 5
BE-4460 Grâce-Hollogne
Belgium
Tel. +32 4 296 82 00
Fax +32 4 296 82 99

ANNEX A for EN 13795-1: Gowns or drapes

Document to join to the samples (1 separate form per reference)

To perform the quoted tests and to provide the correct information in the final report, we would like to receive the following information:

Offer no.

Sample identification

Attention: as this information will be used in the test report, no changes are allowed, once the report is issued.

Tested area(s)/material(s) (if applicable)

Product type (tick your choice)

surgical gown
surgical drape

Identification of sides and manufacturing direction

clearly mark on the sample or describe below the exterior side (= side facing the sterile field):

clearly mark on the sample or describe below the manufacturing direction of the fabric (= weft/cross or chain/machine):



Sterilisation (tick your choice)¹

Note that standard EN 13795-1 requires that all the tests be performed on the sterile finished product (if applicable and except for microbial cleanliness).

In the event the samples do not comply with this requirement, clearly indicate that they are to be tested non-sterile and that you take full responsibility for deviating from it. (*)

- non-sterile sample for microbial cleanliness
- sample is already sterile
- sample can be steam sterilised at 121°C for 15 minutes
- sample can be steam sterilised at 134°C for 4 minutes
- to be tested non-sterile (*)

Precisions on test area

(critical and less critical areas, if necessary, include a picture or drawing)

Important Notice: Samples should be shipped in individual and sealed packages.

Primary contact person for the tests:

Company

Mr./Mrs./Ms.

Phone

Email

Signed for approval

Date:

¹ Some tests require sterilization prior to testing